

MH Hospital
Laboratory Report

bioMérieux Customer:

System #:

Isolate: 01NBC01-1 (Approved)

Card Type: GP Bar Code: 2422669403010083 Testing Instrument: 000015F23819 (11496)

Setup Technologist: Laboratory Supervisor(LabSuper)

Autoprint

Bionumber: 156002365771671

Organism Quantity:

Selected Organism: Enterococcus faecalis

Comments:	

Identification Information	Card: GP	Lot Number: 2422669403	Expires: Mar 18, 2025 12:00 PKT
	Status: Final	Analysis Time: 5.07 hours	Completed: Jun 4, 2022 14:55 PKT
Organism Origin	VITEK 2		
Selected Organism	98% Probability Enterococcus faecalis Bionumber: 156002365771671 Confidence: Excellent identification		
Analysis Organisms and Tests to Separate:			
Analysis Messages:			
Contraindicating Typical Biopattern(s) Enterococcus faecalis LeuA(79),			

Biochemical Details																	
2	AMY	+	4	PIPLC	-	5	dXYL	-	8	ADH1	+	9	BGAL	-	11	AGLU	+
13	APPA	-	14	CDEX	+	15	AspA	+	16	BGAR	-	17	AMAN	-	19	PHOS	-
20	LeuA	-	23	ProA	-	24	BGURr	-	25	AGAL	-	26	PyrA	+	27	BGUR	-
28	AlaA	+	29	TyrA	+	30	dSOR	-	31	URE	-	32	POLYB	+	37	dGAL	+
38	dRIB	+	39	ILATk	-	42	LAC	+	44	NAG	+	45	dMAL	+	46	BACI	+
47	NOVO	+	50	NC6.5	+	52	dMAN	+	53	dMNE	+	54	MBdG	-	56	PUL	-
57	dRAF	-	58	O129R	+	59	SAL	+	60	SAC	+	62	dTRE	+	63	ADH2s	+
64	OPTO	+															

Installed VITEK 2 Systems Version: 9.04.4

MIC Interpretation Guideline:

AES Parameter Set Name:

Therapeutic Interpretation Guideline:

AES Parameter Last Modified:

MH Hospital
Laboratory Report

bioMérieux Customer:

System #:

Isolate: 01NBC04-1 (Approved)

Card Type: GP Bar Code: 2422669403009799 Testing Instrument: 000015F23819 (11496)

Setup Technologist: Laboratory Supervisor(LabSuper)

Autoprint

Bionumber: 010147265763571

Organism Quantity:

Selected Organism: Enterococcus faecium

Comments:	

Identification Information	Card: GP	Lot Number: 2422669403	Expires: Mar 18, 2025 12:00 PKT
	Status: Final	Analysis Time: 4.82 hours	Completed: Jun 4, 2022 14:40 PKT
Organism Origin	VITEK 2		
Selected Organism	87% Probability Enterococcus faecium Bionumber: 010147265763571 Confidence: Acceptable identification		
Analysis Organisms and Tests to Separate:			
Analysis Messages:			
Contraindicating Typical Biopattern(s)			
Enterococcus faecium BGURr(1),AMY(96),BGUR(1),			

Biochemical Details																	
2	AMY	-	4	PIPLC	-	5	dXYL	-	8	ADH1	+	9	BGAL	-	11	AGLU	-
13	APPA	-	14	CDEX	-	15	AspA	-	16	BGAR	+	17	AMAN	-	19	PHOS	-
20	LeuA	-	23	ProA	-	24	BGURr	+	25	AGAL	+	26	PyrA	+	27	BGUR	+
28	AlaA	-	29	TyrA	+	30	dSOR	-	31	URE	-	32	POLYB	+	37	dGAL	+
38	dRIB	+	39	ILATk	-	42	LAC	+	44	NAG	+	45	dMAL	+	46	BACI	+
47	NOVO	(-)	50	NC6.5	+	52	dMAN	+	53	dMNE	+	54	MBdG	+	56	PUL	-
57	dRAF	+	58	O129R	-	59	SAL	+	60	SAC	+	62	dTRE	+	63	ADH2s	+
64	OPTO	+															

Installed VITEK 2 Systems Version: 9.04.4
 MIC Interpretation Guideline:
 AES Parameter Set Name:

Therapeutic Interpretation Guideline:
 AES Parameter Last Modified: