

MH Hospital
Laboratory Report

bioMérieux Customer:

System #:

Isolate: 09NBC01-1 (To be reviewed)

Card Type: GN Bar Code: 2412528503745756 Testing Instrument: 000015F23819 (11496)

Setup Technologist: Laboratory Supervisor(LabSuper)

Printed by: LabSuper

Bionumber: 0007610440104210

Organism Quantity:

Selected Organism: Salmonella ser.Typhi

Comments:	

Identification Information	Card: GN	Lot Number: 2412528503	Expires: Oct 28, 2024 12:00 PKT
	Status: Final	Analysis Time: 5.85 hours	Completed: Jun 15, 2022 20:40 PKT

Organism Origin	VITEK 2
-----------------	---------

Selected Organism	95% Probability Bionumber: 0007610440104210	Salmonella ser.Typhi Confidence: Very good identification
-------------------	------------------------------------------------	--------------------------------------------------------------

Analysis Organisms and Tests to Separate:

Analysis Messages:
 Confirm by serological tests

Contraindicating Typical Biopattern(s)
 Salmonella ser.Typhi PHOS(99),

Biochemical Details																	
2	APPA	-	3	ADO	-	4	PyrA	-	5	IARL	-	7	dCEL	-	9	BGAL	-
10	H2S	-	11	BNAG	-	12	AGLTp	-	13	dGLU	+	14	GGT	+	15	OFF	+
17	BGLU	-	18	dMAL	+	19	dMAN	+	20	dMNE	+	21	BXYL	-	22	BAlap	-
23	ProA	-	26	LIP	-	27	PLE	-	29	TyrA	-	31	URE	-	32	dSOR	+
33	SAC	-	34	dTAG	-	35	dTRE	+	36	CIT	-	37	MNT	-	39	SKG	-
40	ILATk	+	41	AGLU	-	42	SUCT	-	43	NAGA	-	44	AGAL	-	45	PHOS	-
46	GlyA	-	47	ODC	-	48	LDC	+	53	IHISa	-	56	CMT	+	57	BGUR	-
58	O129R	+	59	GGAA	-	61	IMLTa	-	62	ELLM	-	64	ILATa	-			

Installed VITEK 2 Systems Version: 9.04.4

MIC Interpretation Guideline:

AES Parameter Set Name:

Therapeutic Interpretation Guideline:

AES Parameter Last Modified: